

Pressure platforms



sigma
ELECTRONIC PRESSURE PLATFORM



NEW DESIGN



NEW



NEW DESIGN

Platforms	RS	WRS	XL
Active surface	40 x 40 cm	50 x 50 cm	50 x 150 cm
Resistive sensors	1 600 (10 x 10 mm)	4 096 (8 x 8 mm)	12 288 (8 x 8 mm)
Thickness	9 mm	9 mm	9 mm
Built-in WiFi		✓ Batteries included	
USB 2.0 Connection	✓	✓	✓
External dimensions	59 x 56 mm	69 x 61 mm	65 x 161 cm
Weight	3.25 kg	3.8 kg	12 kg
Reference	9760RS	9760WRS	9833XLS
Compatible with	 Podosat 6 , patients software - ref. 9849		



	Ref.	Description
Accessories (sold separately)	9758HL	Protective pouch for SIGMA RS
	9832H	Protective pouch for SIGMA XL
	9750P	Heel positioning system
	9754S	Transparent plexiglas support