

Thermoforming



Remote control and lateral control

Imprinters combined with vacuum press



Foot print



Molding WITHOUT PATIENT

3 times stronger than another system

Highlights

- Molding without patient : **time saver and more comfortable for the patient**
- Perfect thermoforming and perfect molding of the insole
- High temperature thermoforming possibility

FEET PRINT 4.0 EXPERT

Imprinters combined with vacuum press




FEET PRINT 4.0 EXPERT

Ref. 7975B



FEET PRINT 4.0 EXPERT
INTEGRATED COOLING SYSTEM

Ref. 7975BV



Technical characteristics

- 2 feet imprinters
- Vacuum press integrated
- Dimensions : L 915 x W 695 x H 55 mm - Weight : 27 kg

Accessories sold separately



Stand with support handle

Ref. 7975-CP

Dimensions : L 1060 x W 845 mm



Ref. 7613 - **Silicon protection**



Ref. 11110 - **Bactericidal cleaning foam**

Ref. 7909 - **Silicon cleaning spray**