




Thermoforming



Create your own set up



Lateral control
with light indicators



Compensation step
in option

Highlights

- Work from the front or from the back of the patient
- Minimal size
- Easy control and correction

FEET PRINT 3.0

The most ergonomic



FEET PRINT 3.0
1 foot imprinter
Ref. 7979B



Compensation step
PMMA
Ref. 7979MC

Technical characteristics

- Vacuum pump integrated
- Lateral and remote control
- Silicon bag
 - Dual hardness
 - Easy correction
 - Great length useful 420 mm
- Dimensions and weight
 - **Imprinter** : L 540 x W 240 x H 120 mm - 6 kg
 - **Compensation step** PMMA : L 540 x W 200 x H 85 mm - 2 kg

Accessories sold separately



Stand with support handle
Ref. 7979-CP
Dimensions : L 790 x W 790 mm



Ref. 7613 - **Silicon protection**



Ref. 11110 - **Bactericidal cleaning foam**

Ref. 7909 - **Silicon cleaning spray**