

Thermoforming



2.1



2.2

FEET PRINT 2.1 2.2



FEET PRINT 2.1
Ref. 7997B-1B



FEET PRINT 2.2
Ref. 7997B-2B

Technical characteristics

- Vacuum pump integrated
- Remote control
- Silicon bag
 - Dual hardness
 - Easy correction
 - Great length useful 420 mm
- Dimensions : L 690 x W 640 x H 75 mm
- Weight : Feet Print 2.1 : 15 kg - Feet Print 2.2 : 22 kg



2.1 Integrated compensation step



2.2 Independent molding

right foot / left foot
Time saver

Highlights

- Easy control
- Accurate print
- Modern design and easy to clean

Accessories sold separately



Stand with support handle

Ref. 7997-CP

Dimensions : L 900 x W 790 mm



Ref. 7613 - **Silicon protection**



Ref. 11110 - **Bactericidal cleaning foam**

Ref. 7909 - **Silicon cleaning spray**