

Gouge blades



100% made in France
in our workshops in Ecuisses since 1998



500 GOUGES

The dispenser box

Ref. 13002S

Blades n°1 - 2 - 3

Blades n°5 - 8 - 10 - 12



125 GOUGES

The dispenser box

Ref. 13002L

Blades n°1 - 2 - 3

Blades n°5 - 8 - 10 - 12



50 GOUGES

The box

Ref. 13002

Blades n°1 - 2 - 3

Blades n°5 - 8 - 10 - 12



Watch the video of gouge blades manufacturing on our channel.

- ✓ 1 - Cutting of stainless steel, medical grade
- ✓ 2 - Polishing
- ✓ 3 - Sharpening
- ✓ 4 - Disinfection and packaging in a clean room
- ✓ 5 - Sterilization



Ref.	Description	
13002M	50 foam gouge blades Steriles,posables	<i>The box</i>
13000	Gouge blade Autoclavable, removable and perforated	<i>Per unit</i>