

Unit

Choose your front color



Example : right hand
Tray : black
Light : SUNLED
Front color : orange U340



Example : left hand
Tray : white
Light : MAIA PODO
Front color : grey U963



Tools support

- Integrated to the workplan
- Automatic detection of tools
- Cleared workplan



Workplan

- Laminate surface
- Glass protection, white or black



Specific drawers

- High quality runners with shock absorbers
- Removable interiors fittings




Easy move

- Castors with brakes
- Easy maintenance of grounds



Footcontrol

- Easy use
- Regulation of the micromotor and turbine speed
- Removable footcontrol



Connection to water supply

- **Free (replace the tank)**
- Possibility to plug the unit straight to the water network
- Demineralizing cartridge
- No stagnant water, hygienic - No tank filling

NEW



HARMONIE

Ref. 10850H

Configuration

- Workplan - white or black glass
- Tool kit not included - check page 124
- Light not included - check page 140
- Compressor not included - check page 141