

Unit



Example : right hand
Tray : white glass
Light : SUNLED
Arm system : H-Flex - 3 HC



H-Flex by Capron

- Great length
- Adapted for left and right hand
- Automatic detection of tools
- Available with 2, 3 or 4 tools



Springflex by Capron

- Soft manipulation in all directions
- Easy maintenance of the glass protection
- Automatic detection of tools
- Available with 2, 3 or 4 tools



CLOUD

Ref. 10810C

Configurations

- Tool kit not included - check page 124
- Light not included - check page 140
- Compressor not included - check page 141

In option

Workplan

Ceramic or Quartz



Example : right hand with suction
Tray : white glass
Light : MAIA PODO
Arm system : Springflex 4 tools



Swivel Arm

- Easy manipulation
- Adjustable arm control




Compensated arm control

- Articulated and compensated arm on 4 axis
- Height adjustment
- Facilitate the work clearance
- Bring the instruments really close to the foot (right foot - left foot)



Foot control

- Progressive regulation of micromotor and turbine
- Separated spray control
- Removable foot control



Connection to water supply

- **Free (replace the tank)**
- Possibility to plug the unit straight to the water network
- Demineralizing cartridge
- No stagnant water, hygienic - No tank filling