

Chairs

Customizable
covering colors



1 COLOR
20287



1 COLOR
20280




BI COLOR 1
00002 - 20299



BI COLOR 2
00002 - 20299



Sliding leg rests
with double extension : +33 cm



Reclining leg rests
for better comfort

EXCLUSIVITY CAPRON

Our concept of tilting legrest **relaxes the patient's leg** and foot.
Results : More comfortable for the patient and for the doctor.



Laterally fold back leg rests
Front access to the patient



Foot control
2 memories and automatic reset

NEW



GRAPHITE

Ref. 11670

Characteristics

3 independent movements	
Motors.....	2
Up and down.....	electric
Trendelenburg.....	electric
Back.....	manual - gas lift
Rotation.....	360° with locking
Low position.....	56 cm
High position at leg end.....	139 cm
Width of seat.....	52 cm
Lifting capacity.....	220 kg
Weight.....	110 kg



Accessories sold separately

Ref.	Description
11650T	Removable headrest
10950-A	Disposable dust tray support
10966	Plastic dust tray
10964	2 protections for leg rests

NEW