

# Tailor-Fit

(MODEL 49)


**3**  
YEAR  
GUARANTEE

 Customised  
Just For YOU

Vest with 15cm partial overlap, with choice of zipper or Velcro fastening.



Velcro closure



Zipper closure

## PRODUCT OVERVIEW

- ✓ Greater air circulation
- ✓ Padded shoulders
- ✓ Soft outer material to allow for complete flexibility
- ✓ Even transfer of weight between shoulders and hips to reduce fatigue
- ✓ Made to individual measurements
- ✓ Antibacterial surface

## AmRay OPTIONS

- ✓ Monogram
- ✓ Logo on pocket
- ✓ Thyroid collar clipped on



## MATERIALS

Light Lead

Light Plus

Supreme Lead-Free

## LEAD EQUIVALENCE

.25/.25

.35/.25

.35/.35

.50/.25

.50/.35

## PROTECTION



Front, sides and  
back protection  
Partial front overlapping

## SIZE OPTIONS

XS

S

M

L

XL

2XL

Other sizes available on request  
Please refer to size chart p. 12