

Ergo-Fit+

(MODEL 45)



3
YEAR
GUARANTEE



Customised
Just For YOU

Full protection to front and back with adjustable elasticated Velcro flaps for back support and individual comfort.



Velcro flaps

PRODUCT OVERVIEW

- ✓ Greater air circulation
- ✓ Padded shoulders
- ✓ Soft outer material to allow for complete flexibility
- ✓ Even transfer of weight between shoulders and hips to reduce fatigue
- ✓ Made to individual measurements
- ✓ Antibacterial surface

AmRay OPTIONS

- ✓ Monogram
- ✓ Logo on pocket
- ✓ Thyroid collar clipped on
- ✓ Velcro or buckle closure



MATERIALS

Light Lead

Light Plus

Supreme Lead-Free

LEAD EQUIVALENCE

.25/.25

.35/.25 .35/.35

.50/.25 .50/.35

PROTECTION



Front, back
and sides protection

SIZE OPTIONS

XS S M

L XL 2XL

Other sizes available on request
Please refer to size chart p. 12