

Easy-Fit

(MODEL 44)



Double sided vest with Velcro closure to front.



Velcro closure – Shoulder



Velcro closure – Side

PRODUCT OVERVIEW

- ✓ Greater air circulation
- ✓ Padded shoulders
- ✓ Soft outer material to allow for complete flexibility
- ✓ Even transfer of weight between shoulders and hips to reduce fatigue
- ✓ Made to individual measurements
- ✓ Antibacterial surface

AmRay OPTIONS

- ✓ Monogram
- ✓ Logo on pocket
- ✓ Thyroid collar clipped on



MATERIALS

Light Lead

Light Plus

Supreme Lead-Free


LEAD EQUIVALENCE

.25/.25 .35/.25

.35/.35 .50/.25

.50/.35 .50/.50

PROTECTION



Front and back protection

SIZE OPTIONS

XS

S

M

L

XL

2XL

Other sizes available on request
Please refer to size chart p. 12