


Reinforced Silicone Laryngeal Mask Airway 加强型硅胶喉罩



1. Made of imported medical grade silicone, high quality.
2. Silicone cuff is made by one time, provides good and soft sealing.
3. Flexible tube with spiral reinforced wire, reducing the risk of kinking.
4. Single-use and reusable 2 types available.
5. Reusable type is autoclave sterilization only, can be reused up to 40 times.

Size	1.0	1.5	2.0	2.5	3.0	4.0	5.0
Single use Item No.	LM1004S1R	LM1507S1R	LM2010S1R	LM2514S1R	LM3020S1R	LM4030S1R	LM5040S1R
Single use Item No. (With Bar)	LM1004GS1R	LM1507GS1R	LM2010GS1R	LM2514GS1R	LM3020GS1R	LM4030GS1R	LM5040GS1R
Reusable Item No.	LM1004S2R	LM1507S2R	LM2010S2R	LM2514S2R	LM3020S2R	LM4030S2R	LM5040S2R
Reusable Item No. (With Bar)	LM1004GS2R	LM1507GS2R	LM2010GS2R	LM2514GS2R	LM3020GS2R	LM4030GS2R	LM5040GS2R
Patient Selection	0-5KG	5-10KG	10-20KG	20-30KG	30-50KG	50-70KG	70-100KG
Max. Inflation Volume	4ml	7ml	10ml	14ml	20ml	30ml	40ml
Remark	1. Bar/Aperture: With <input checked="" type="checkbox"/>		Without <input checked="" type="checkbox"/>				
	2. Cuff Color: Flesh <input checked="" type="checkbox"/>		Blue <input checked="" type="checkbox"/>		Transparent <input checked="" type="checkbox"/>		

