



## Skirt

(MODEL 04)



MATERIALS	LEAD EQUIVALENCE	SIZE OPTIONS	PROTECTION
Light Lead	.35/.25 .35/.35	Size S M L	<p>Front, back and sides protection</p>
Light Plus	.50/.25 .50/.35	Age 0-4 yrs 5-8 yrs 9-12 yrs	
Supreme Lead-Free	.70/.25 .70/.35	Width 50 54 59	
		Length 35 40 50	

Note: Child motif subject to change

## Dental Apron /OPG

(MODEL 10)

### PRODUCT OVERVIEW

- ✓ Practical Velcro fastening for ease of fitting and removal

MATERIALS	LEAD EQUIVALENCE	PROTECTION	SIZE OPTIONS
Light Lead	.25mm Pb (Only available for Light)	<p>Front protection</p>	Width: 40cm
Light Plus	.55mm Pb		Length: 70cm
Supreme Lead-Free	.5mm Pb		Other sizes available on request

## Ovarian Shield

(Child)



PRODUCT OVERVIEW	MATERIALS
<ul style="list-style-type: none"> <li>✓ Paediatric Ovarian Shield</li> <li>✓ Set of 4</li> <li>✓ Velcro Belt for attachment</li> <li>✓ Sealed edges for hygiene control</li> </ul>	Light Lead
SIZE AVAILABLE	Light Plus
S M L	Supreme Lead-Free
LEAD EQUIVALENCE	
1mm Pb	