

Paediatric

Effective radiation protection is achieved when garments are made to the correct size. **AmRay** offers the widest choice of protective aprons for children. A full choice of colour options are available as shown on our colour swatch.



Coat Apron

(MODEL 06)

MATERIALS	LEAD EQUIVALENCE	SIZE OPTIONS	PROTECTION																
Light Lead	.25/.25	<table border="1"> <tr> <th>Size</th> <th>S</th> <th>M</th> <th>L</th> </tr> <tr> <th>Age</th> <td>0-4 yrs</td> <td>5-8 yrs</td> <td>9-12 yrs</td> </tr> <tr> <th>Width</th> <td>42-48</td> <td>44-50</td> <td>47-55</td> </tr> <tr> <th>Length</th> <td>55</td> <td>60</td> <td>70</td> </tr> </table>	Size	S	M	L	Age	0-4 yrs	5-8 yrs	9-12 yrs	Width	42-48	44-50	47-55	Length	55	60	70	<p>Front, back and sides protection</p>
Size	S	M	L																
Age	0-4 yrs	5-8 yrs	9-12 yrs																
Width	42-48	44-50	47-55																
Length	55	60	70																
Light Plus	.35/.25 .35/.35																		
Supreme Lead-Free	.50/.25 .50/.35																		
<p>Note: Child motif subject to change</p>																			



Double-Sided Apron

(MODEL 05)

MATERIALS	LEAD EQUIVALENCE	SIZE OPTIONS	PROTECTION
Light Lead	.25/.25	<p>Width: 40cm</p> <p>Length: 67cm</p> <p>Other sizes available on request</p>	<p>Front and back protection</p>
Light Plus	.35/.25 .35/.35		
Supreme Lead-Free	.50/.25 .50/.35		



Single Front Apron

(MODEL 07)

- Closure A:** Velcro flap between the shoulders
- Closure B:** Buckle with fully adjustable strap
- Closure C:** Elasticated Velcro fastening provides individual fit

MATERIALS	LEAD EQUIVALENCE	SIZE OPTIONS	PROTECTION
Light Lead	.25mm Pb <small>(Only available for Light)</small>	<p>Width: 40cm</p> <p>Length: 67cm</p> <p>Other sizes available on request</p>	<p>Front protection</p>
Light Plus	.35mm Pb		
Supreme Lead-Free	.5mm Pb		