

Patient Protection

There has been a significant increase in the use of ionizing radiation in medicine during recent decades, resulting in patients being increasingly exposed. Our comprehensive range of accessories ensures maximum protection and comfort for the patients.



Triangle



PRODUCT OVERVIEW		MATERIALS
✓ Male Gonad Protector		Light Lead
SIZE AVAILABLE	LEAD EQUIVALENCE	Light Plus
22cm sides	1mm Pb	Supreme Lead-Free
✓ Sealed edges for hygiene control		

Gonad Shield



PRODUCT OVERVIEW		MATERIALS
✓ Male Gonad Protector ✓ Sealed edges for hygiene control		Light Lead
SIZE AVAILABLE	LEAD EQUIVALENCE	Light Plus
S: 0-9 yrs	1mm Pb	Supreme Lead-Free
M: 9-14 yrs		
L: 14+		

Ovarian Shield (Adult & Child)



PRODUCT OVERVIEW		MATERIALS
✓ Female Ovarian Shield ✓ Set of 4 ✓ Velcro Belt for attachment ✓ Sealed edges for hygiene control		Light Lead
SIZE AVAILABLE	LEAD EQUIVALENCE	Light Plus
S M L	1mm Pb	Supreme Lead-Free
Also available in paediatric size.		