



Blanket



PRODUCT OVERVIEW

- ✓ Versatile and flexible patient protection

MATERIALS

Light Lead

Light Plus

Supreme Lead-Free

LEAD EQUIVALENCE

0.25mm Pb

0.35mm Pb

.5mm Pb

SIZES

Sizes available on request

Portable Shield

PRODUCT OVERVIEW

- ✓ The Adjustable Mobile Shield can be set between 91cm and 152 high
- ✓ The nylon covered .50mm Pb shield is easy to maneuver with its sturdy, 26" star-leg base and durable castor wheels

