

OVERBED TABLE - ADJUSTABLE HEIGHT

MODEL WP2100



Product Features

- Steel construction
- Scratch & chip resistant epoxy powder coating
- Height adjustable classic tension spring Knob lock
- Laminated Chipboard surface with all round T-Frame rubber trimmings around table top edge
- Mobile on 4 x 50mm dia. swivel castors

Technical Specifications

- Table Top Approx. Dim(mm): 813(L) x 356(W)
- Overall Approx. Dim(mm): 825(L) x 406(W) x 760-1080(H)