




Double sided drawers make supplying medicines while on the go fast and efficient.

- PP worktop and handles
- Aluminium frame
- Central key lock control on both sides
- ABS side panels
- 4" TPR casters (2 w/ brake)
- 3" bin drawer x 2 rows per side
3" ABS bin x 20pcs
(5pcs per row x 2 rows x 2 sides)
- 6" bin drawer x 1 row per side
6" ABS bin x 10pcs
(5pcs per row x 1 row x 2 sides)
- 6" drawer x 1 + 9" drawer x 1



FC3432DXL

ABS Bin



ABS bins with PS dividers come in 2 sizes and easily separate medicines and supplies for just about any need.

- 1 3" Height (W 112 x D 215 x H 75 mm)
- 2 6" Height (W 112 x D 215 x H 154 mm)



Mount accessories on both sides.

SPECIFICATIONS

FC3432DXL

	mm	inch
W	760	29.9
D	490	19.3
H	1090	42.9
H (w/o 4" caster)	990	39.1



FC3440DXL



ABS bins are easy to remove and replace.

Drawers



Strong drawers with large spaces easily accommodate large expendables and supplies.

- PP worktop and handles
- Aluminium frame
- Central key lock control on both sides
- ABS side panels
- 4" TPR casters (2 w/ brake)
3" bin drawer x 2 rows per side
3" ABS bin x 20pcs
(5pcs per row x 2 rows x 2 sides)
- 6" bin drawer x 2 rows per side
6" ABS bin x 20pcs
(5pcs per row x 2 rows x 2 sides)

SPECIFICATIONS

FC3440DXL

	mm	inch
W	760	29.9
D	490	19.3
H	1090	42.9
H (w/o 4" caster)	990	39.1