

## DOUBLE MEDICAL CART

This double medical cart provides large space for storage and organizes medicines easily and clearly. The key lock design secures all drawers. It also offers stable, smooth mobility and enhanced convenience with changeable accessories to meet your needs.

### VTB-2-DX

(18 tilt bins (S)+1010.5mm rail(VPR-4)x2)



### VBS-DX

(Bulk storage bin x 9+one metal hanger +1010.5mm rail (VPR-4)x1)



(Door w/shelves is available)



DX34

### A drawer dimension

- W569 x D398 x H76mm (3")
- W569 x D398 x H155mm (6")
- W569 x D398 x H234mm (9")
- 3" divider sets: P-FD3 / P-AD3
- 6" divider sets: P-FD6 / P-AD6

### B drawer dimension

- W399 x D398 x H76mm (3")
- W399 x D398 x H155mm (6")
- W399 x D398 x H234mm (9")
- 3" divider sets: N-FD3 / N-AD3
- 6" divider sets: N-FD6 / N-AD6

### Drawer combination

You can select different drawer combinations within the limited height 27".

### SPECIFICATIONS

DX34

	mm	inch
W	1125	44.3
D	545	21.5
H	1010	39.8
H (w/o 5" Casters)	885	34.8

### Features

- Metal construction with powder coating
- 4 aluminum frames
- ABS top, pull-out work table & push handle
- Central key lock design
- Standard 5" iron casters  
(1 free, 1 w/direction control, 2 w/brake)



VLP

VVP

VLN

### SCREEN HOLDER


If using PC monitor, you can choose VVP or VLP to assemble on the right side of carts. It is height adjustable.

- Part no. VLP includes 2 + 3
- Part no. VVP includes 1 + 2

### NOTEBOOK HOLDER

If using the notebook, you can choose

- Part no. VLN includes 2 + 3 + 4
- Part no. VVN includes 1 + 2 + 4



### 1 VESA connector

You can assemble PC screen or notebook holder on this connector to make your work efficient.

- Weight Capacity: ~8kg (17.6 lbs)
- Rotation: 270°
- Pivot: 180°
- VESA: 75/100mm

### 2 Support holder

70cm length, metal with powder coating, height adjustable by assembling in different positions

### 3 Monitor long arm

This monitor long arm provides you the flexible distance when working.

- Max weight: 10kg,
- Max length: 25cm,
- VESA: 75/100mm

### 4 Notebook holder

Height adjustable