


## ✚ Standard Endotracheal Tube (Cuffed)



1. Made of medical grade PVC, clear and smooth.
2. High volume, low pressure cuff maintains good sealing.
3. Murphy eye design, tube with X-ray opaque line.

Item No.	Size ID(mm)	Item No.	Size ID(mm)
ET30PC	3.0	ET65PC	6.5
ET35PC	3.5	ET70PC	7.0
ET40PC	4.0	ET75PC	7.5
ET45PC	4.5	ET80PC	8.0
ET50PC	5.0	ET85PC	8.5
ET55PC	5.5	ET90PC	9.0
ET60PC	6.0	ET95PC	9.5

Notes: Uncuffed type is available, from ID2.0-ID9.5

## ✚ Reinforced Endotracheal Tube (Cuffed)

1. Made of medical grade PVC, clear and smooth.
2. High volume, low pressure cuff maintains good sealing.
3. Tube with spiral reinforced wire reducing risk of kinking.
4. Murphy eye design, tube with X-ray opaque line.

Item No.	Size ID(mm)	Item No.	Size ID(mm)
RET30PC	3.0	RET65PC	6.5
RET35PC	3.5	RET70PC	7.0
RET40PC	4.0	RET75PC	7.5
RET45PC	4.5	RET80PC	8.0
RET50PC	5.0	RET85PC	8.5
RET55PC	5.5	RET90PC	9.0
RET60PC	6.0	RET95PC	9.5



Notes: can be equipped with stylet. Uncuffed type is from ID2.0-ID9.5

