

ISOLATION CARTS

MX Isolation Carts were designed with bright yellow drawers to alert staff and visitors. This series is available with non-locking drawers, central locking drawers, or individually locking drawers. Typical accessory options include deep drawers for bulky supplies (gowns, caps, gloves, etc.) and a built-in aluminium frame with yellow sign line.



SPECIFICATIONS	MX28ISO		MX31ISO		MX34ISO		MX37ISO	
	mm	inch	mm	inch	mm	inch	mm	inch
W	690	27.2	690	27.2	690	27.2	690	27.2
D	518	20.4	518	20.4	518	20.4	518	20.4
H	848	33.4	927	36.5	1006	39.6	1085	42.7
H (w/o 5" Casters)	714	28.1	793	31.2	872	34.3	951	37.4



Drawers & dividers
Strong drawers with large storage capacity easily store large expendables and other supplies.



Accessory combinations
Height adjustable accessories can be set up to suit any requirement.

STANDARD ACCESSORIES

- 1 Multi storage (L)
- 2 Waste container w/ lid
- 3 Glove dispenser



Multi storage (L)



Waste container w/ lid



Glove dispenser