

LAB-X SERIES, BLOOD DRAWING CHAIRS

Clinton's Lab X Series, Blood Drawing Chairs put patient comfort first. A thick, padded, contoured seat and arched backrest are supported by a stable, durable, steel, powder-coated frame. Height-adjustable arms make operation easy for clinician and comfortable for patient.

Lab X Series Chairs feature:

- 11/4" square, heavy duty, all-welded, tubular steel, frame (3 cm)
- Adjustable height and depth flip arm
- Flip arm positions in place when in use, then flips out of the way for easy exit
- Optional ClintonClean™ Arms available
- Comfortable, 2" padding throughout (5 cm)
- Neutral Gray, powder-coated frame
- Height adjustable armrest
- Easy-clean, durable vinyl in a selection of colors
- Chair-mounted options available (see page 21)
- 400 lbs. load capacity under normal use (181.4 kg)



Armrests are adjusted by two star knobs



Adjustable depth flip-arm



66010

Lab X Series, Blood Drawing Chair with Padded Arms

- Contoured padded seat and padded, arched backrest
- Padded flip arm and armrests
- See page 21 for chair-mounted options
- Option -665, ClintonClean™ solid plastic arms

66010

dimensions		seat height	
seat	17 1/2" x 18"	19"	48 cm
overall	width 31"	depth 26"	height 36"
	79 cm	66 cm	91 cm
arms	arms height range 6" 15 cm		

IMPORTANT NOTE: Clinton Blood Drawing Chairs are not engineered to restrain a fainting or falling occupant. The chair's occupant should NEVER be left unattended.



-665 ClintonClean™ Arms option

Lab X Series, Blood Drawing Chair with Padded Flip Arm and Drawer

- Contoured padded seat and padded, arched backrest
- Padded flip arm and armrest
- ClintonClean™ drawer top
- One-piece, all-plastic drawer
- See page 21 for chair-mounted options
- Option -665, ClintonClean™ solid plastic arms

66020

dimensions		seat height	
seat	17 1/2" x 18"	20"	50.8 cm
overall	width 38"	depth 27"	height 36"
	97 cm	69 cm	91 cm
arms	arms height range 6" 15 cm		

Lab X Series, Extra-Tall, Blood Drawing Chair with Padded Arms

- Contoured padded seat and padded, arched backrest
- Padded flip arm and armrests
- Extra height reduces bending for phlebotomist
- See page 21 for chair-mounted options
- Option -665, ClintonClean™ solid plastic arms

66011

dimensions		seat height	
seat	17 1/2" x 18"	26"	66 cm
overall	width 31"	depth 27"	height 42 1/2"
	79 cm	69 cm	108 cm
arms	arms height range 6" 15 cm		

Lab X Series, Extra-Tall, Blood Drawing Chair with Padded Flip Arm and Drawer

- Contoured padded seat and padded, arched backrest
- Padded flip arm and armrest
- ClintonClean™ drawer top
- One-piece, all-plastic drawer
- Extra height reduces bending for phlebotomist
- See page 21 for chair-mounted options
- Option -665, ClintonClean™ solid plastic arms

66022

dimensions		seat height	
seat	17 1/2" x 18"	26"	66 cm
overall	width 38"	depth 27"	height 42 1/2"
	97 cm	69 cm	108 cm
arms	arms height range 6" 15 cm		



66020



66011



66022

All blood drawing products are non-returnable unless they are unused in unopened, original packaging.