


PHLEBOTOMY CARTS

Easy to handle and maneuver in tight spaces with ample space for supplies and equipment, Clinton's phlebotomy carts are the right fit for busy nurses and phlebotomists.

5-Leg, Phlebotomy Carts feature:

- Gray laminate cabinet construction
- Easy-clean, vacuum-formed, gray, ABS, plastic top with textured finish and raised rim on all sides
- Lockable door with soft-close hinges
- Removable plastic storage tray(s) on easy slide tracks
- Easy grip ClintonClean™, solid plastic handle
- Pneumatic height adjustment
- Stable, 5-leg, aluminum base
- Smooth rolling, 3" ball bearing casters (8 cm)
- Optional glove box holder available



Optional glove box holder mounts to back (Order Single White Steel Glove Box Holder GW-2000)

One or two plastic bins slide on plastic channels



One-Bin Phlebotomy Cart

67100			
top	width 20" 51 cm	depth 20 ³ / ₄ " 53 cm	
cabinet	width 13 ¹ / ₄ " 34 cm	depth 17 ¹ / ₄ " 44 cm	height 7 ³ / ₄ " 20 cm
	base diameter 24" 61 cm		
height range	33" – 43" 84 cm – 109 cm		

Two-Bin Phlebotomy Cart

67200			
top	width 20" 51 cm	depth 20 ³ / ₄ " 53 cm	
cabinet	width 13 ¹ / ₄ " 34 cm	depth 18" 46 cm	height 10 ¹ / ₂ " 27 cm
	base diameter 24" 61 cm		
height range	33" – 43" 84 cm – 109 cm		



(Tray on top and supplies shown not included)

H-Base, Phlebotomy Carts feature:

- Gray laminate cabinet construction
- Easy-clean, vacuum-formed, gray, ABS, plastic top with textured finish and raised rim on all sides
- Lockable door with soft-close hinges
- Removable plastic storage tray(s) on easy-slide tracks
- Easy grip ClintonClean™, solid plastic handle
- Pneumatic height adjustment
- H-leg, powder-coated steel base
- Smooth rolling, 3" ballbearing casters (8 cm)
- Optional glove box holder available

Small, H-Base, Pneumatic, Two-Bin Phlebotomy Cart

67202			
top	width 20" 50.8 cm	depth 20 ¹ / ₂ " 52 cm	
cabinet	width 13 ¹ / ₄ " 33.65 cm	depth 18" 45.72 cm	height 10 ¹ / ₂ " 26.67 cm
	base diameter 24" 61 cm		
height range	29" – 42 ³ / ₄ " 73.66 cm – 108.58 cm		



67202

Two plastic bins slide on plastic channels



1" powder-coated, steel tubing handle

(Tray on top and supplies shown not included)

Large, H-Base, Pneumatic, Four-Bin Phlebotomy Cart

67203			
top	width 20" 51 cm	depth 20 ¹ / ₂ " 52 cm	
cabinet	width 13 ¹ / ₄ " 34 cm	depth 17 ³ / ₄ " 45 cm	height 19 ¹ / ₂ " 50 cm
	base diameter 24" 61 cm		
height range	29" – 42 ³ / ₄ " 74 cm – 109 cm		



67203



Four plastic bins slide on plastic channels

(Tray on top and supplies shown not included)