



# Mobile Protective Screens

Mobile Protective Screens are used to provide optimum protection to medical personnel during procedures. As they are used within various healthcare departments the requirements for a cleanable surface is of the upmost importance. All **AmRay** Protective Screens are constructed and finished with a encapsulated, antibacterial surface which is both durable and cleanable.

## AMS3040

### PRODUCT OVERVIEW

- ✓ Multipurpose mobile screen
- ✓ Completely seamless hygienic, easily cleanable PVC finish
- ✓ Easy move castors with brake
- ✓ 2mm lead core supplied as standard (alternative thickness available)
- ✓ Unique Flush Vision panel design with 2mm Pb lead equivalence
- ✓ Standard size 83cmx190cm (WxH)
- ✓ Standard glass size 30x40cm
- ✓ Other sizes available on request



## AMS6868

### PRODUCT OVERVIEW

- ✓ Multipurpose mobile screen
- ✓ Completely seamless hygienic, easily cleanable PVC finish
- ✓ Easy move castors with brake
- ✓ 2mm lead core supplied as standard (alternative thickness available)
- ✓ Unique Flush Vision panel design with 2mmPb lead equivalence
- ✓ Standard size 83cmx190cm (WxH)
- ✓ Standard glass size 68x68cm
- ✓ Other sizes available on request