

## **HOSPITAL ELECTRICAL RECLINING CHAIR**

### **MODEL CCH0042**



#### **Product Features**

- Upholstery with PVC cover with PU foam
- Reclining seat and leg rest operated via a remote control
- Mobile on 4 x 50mm dia. castors with 2 brakes

#### **Technical Specifications**

- Backrest Height Approx. Dim (mm): 635mm
- Seat Approx. Dim (mm): 508(W) x 559(D) x 483(H)
- Sleeping area Approx. Dim (mm): 508(W) x 1778(D)
- Approx. Dim (mm): 762(W) x 800(D) x 978(H)