

RECLINING CHAIR **MODEL CCH0039-3**



Product Features

- Steel frame construction
- Scratch and chip resistant epoxy powder coating
- Upholstery with PVC cover with PU foam
- Backrest adjustable by gas spring and footrest adjustable by friction lock
- Mobile on 4 x 125mm castors with 2 brakes fitted

Technical Specification

- Seat Area Approx. Dim (mm): 540(W) x 465(D) x 470(H)
- Sleeping Area Approx. Dim (mm): 1550(D) x 540(W) x 470(H)
- Backrest Height Approx. Dim (mm): 710
- Approx. Dim (mm): 625(W) x 850(D) x 1100(H)