

RECLINING CHAIR-ELECTRIC

MODEL CCH0039-E



Product Features

- Steel frame work with epoxy powder coated
- Cushion padded with Mediflex cover, fire retardant, anti-bacteria, anti-static and water proof for the ease to cleaning
- Reclining chair by electric motorized
- Backrest and legrest able to adjust separately
- Armrest with cushion padded with Mediflex cover height adjustable and swing able.
- Mobile on 4 x 125mm diameter castors with 2 brakes
- Able to reclining into flat position

Technical Specifications

- Backrest height: 802mm
- Armrest: 178mm(W)
- Overall Approx. Dim(mm): 965(W) x 914(D) x 1245(H)
 - Seat area from floor Approx. Dim (mm): 550(W) x 508(D) x 495(H)
 - Sleeping Area Approx. Dim: 550(W) x 1753(D)
- Reclining position: 180°
- Optional: a) IV pole bracket with IV pole
b) Battery back-up