

ELECTRIC DELIVERY BED
MODEL DB4000



**Spread Footrest
Position**



Lithothomy Position



Product Features

- Steel construction finishing with scratch and chip resistant of hard epoxy powder coating
- Convertible from bed to delivery table
- Hi-Lo, backrest and trendelenburg position controlled by electrical motorized
- Foot controls unit to allow a hand free for medical personnel
- Built in backup battery
- Retractable leg section
- Detachable PE head and foot boards
- Double easy fold-away PE blow side rails with gas spring assisted adjustment
- Manual Override CPR release lever at both sides for backrest
- Integrated with one pair handgrips at both sides of the bed
- Delivery position, Squatting position, spreaded footrest position & lithothomy position
- Removable mattress support base
- Complete with accessories kit stand
- Complete with 3-section PU foam mattress with mediflex cover (water proof, fire retardant, anti-bacterial, anti-static)
- Mobile on 4 x 150mm dia. Tente (Germany) central locking system

Electrical System

- LINAK DC-actuator system (TUV & UL approved & EC60601-1)
- 24VDC out-put
- Built in back-up battery

Technical Specification

- Backrest: 0° - 75°
- Trendelenburg: 0° - 8°
- Overall Approx. Dim (mm): 1494(L) stowed/2137(L) extended x 990(W) x 546-930(H)
- Safety load capacity: 250kg UDL

Accessories Included

- CPR Quick Release Level for Backrest
- Patient Control Handset
- Electrical Foot Control Pedal

Accessories Included

- IV drip pole x 1 unit
- Swivel padded knee crutch x 1 pair
- Spreaded footrest x 1 pair
- Patient step x 1 unit
- Squatting support x 1 pair
- Padded mattress for leg section
(150mm Thick) x 1 unit
- Stainless steel fluid collection bowl x 1 unit
- Accessories kit stand mobile on
75mm diameter castors x 1 unit

Convenient & Space Efficient Mobile Accessory Kit Stand

