

HOSPITAL ELECTRICAL BED **MODEL BA7001**



Product Features

- Ultra-low with the minimum height at 240mm (9.5")
- All parts are made by tooling to ensure uniform quality
- Steel construction finishing with scratch and chip resistant of hard epoxy powder coating
- Electric motorize bed with smooth and control operation by 3 actuators
- Detachable head and foot boards
- Mobile on 4 x 100mm diameter castors with 2 brakes diagonally fitted
- Hi-Lo, backrest & kneerest controlled by hand set control
- 4 provisions for IV pole
- 2 drainage bag holders
- 4 provisions for monkey pole
- With 4 corners rotating bumpers
- Intention of use: suitable to be use in ward and geriatric ward

Technical Specifications

- Hi-Lo: 240-650mm
- Working capacity: 180kg
- Backrest: 0° - 75°
- Auto Contour
- Kneerest: 0° - 25°
- Legrest: 0° - 45°
- Overall Approx. Dim(mm): 2115(L) x 1013(W)

Electrical System

- Linak DC-actuator system (TUV & UL approved & EC60601-1)
- 24 DC Output
- Voltage: 220-240; 50Hz - 60Hz
- Built in backup battery

Accessories Included

- Medi-flex mattress with PU foam (32 Density) x 1 unit (**BX411**)
- Mattress Approx. Dim (mm):1905 (L) x 864 (W) x 102 (H)
- Mattress cover with fire retardant, anti-bacterial, anti-static and waterproof x 1 unit
- Collapsible aluminium side rail x 1 pair; Approx. Dim (mm): 410 (H) to 1440 (L) (**BX200**)

Optional Accessories

- Monkey pulley (**BX910**)
- IV pole with 2 hooks- chrome plated (**BX300**)
- IV pole with 2 hooks- stainless steel (**BX301**)
- Castors with Central Locking System
- Orthopedic Frame