

STEP-UP TREATMENT TABLES

Step-Up, Treatment Table with Doors & Drawers

- Sturdy, all laminate and steel construction
- Handy, steel stool with slip-resistant top pulls out for use, and neatly stores out of the way when not in use
- Stool can be easily changed on site for use on either side of the table
- Adjustable backrest with 2" for firm foam padding
- 3-door storage accessible from either side of the table
- 2-drawer module can be easily moved for access from either left or right side of the table
- Non-marking glides for easy positioning

9705-27	length	width	height
	72"	27"	31"
	183 cm	69 cm	79 cm

9705-30 30" width model (76 cm)

Flat Top Option Use -088 after model number



9705



Flat Top



Stool Stored



9704



Flat Top



Stool Stored


Step-Up, Treatment Table with 4 Doors

- Sturdy, all laminate and steel construction
- Handy, steel stool with slip-resistant top pulls out for use, and neatly stores out of the way when not in use
- Stool can be easily changed on site for use on either side of the table
- Adjustable backrest with 2" for firm foam padding
- 4-door storage accessible from either side of the table
- Non-marking glides for easy positioning

9704-27	length	width	height
	72"	27"	31"
	183 cm	69 cm	79 cm

9704-30 30" width model (76 cm)

Flat Top Option Use -088 after model number

 **WARNING:** Products featured on this page can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer.

For more information go to www.P65Warnings.ca.gov.

