



NEEDLE GUIDE ENDOCAVITY



A vital physician's tool to be used as an accessory to ultrasound transducers, guiding a needle or catheter into a target to perform a biopsy for precise needle placement.

Fits: Esaote EC123, TRT12, IVT12, Philips E6509, GE IC5-9-RS, Samsung Medison BCC4-9, Pie Medical C7-4R8TRD, Siemens BE9-4, Terason 8EC4

Emerald Sterile Disposable

- o One-time use saves time & resources from cleaning and sterilizing reusables
- o Wide set funnel & clear design allows for precise needle placement
- o Smooth profile ensures ultimate patient comfort
- o Attaches securely to probe with ease



Order Information

Works with 16 gauge or thinner

Ref #5145 - Sterile Disposable Needle Guides, 24/box

Contains: Needle Guide Only

Ref #5147 - Sterile Disposable Procedure Kit (latex), 24/box

Contains: Needle Guide, Latex Probe Cover, Sterile Gel, and Bands

Ref #5148 - Sterile Disposable Procedure Kit (not made with latex), 24/box

Contains: Needle Guide, Probe Cover, Sterile Gel, and Bands



For Nearest Distributor & Additional Order Information Contact:

4125 Westcor Court, Coralville, IA 52241 Phone: (319)545-7100 Fax: (319)545-7099

Email: sales@protekmedical.com Web: www.protekmedical.com