

Lead Caps

| PRODUCT OVERVIEW | MATERIALS | LEAD EQUIVALENCE |
|--|-------------------|--|
| <ul style="list-style-type: none"> ✓ The lead cap offers protection for your head against scatter radiation ✓ Size adjustable with elastic closure at the back | Light Lead | <div>.35mm Pb</div> <div>.5mm Pb</div> |
| | Light Plus | |
| | Supreme Lead-Free | |



Customised
Just For YOU

Mesh Cap



Elastic Hat




Lead Gloves

| PRODUCT OVERVIEW | SIZE OPTIONS | LEAD EQUIVALENCE |
|---|--------------|--------------------|
| <ul style="list-style-type: none"> ✓ One-Piece flexible lead vinyl for easy grip ✓ Inner cotton glove supplied for greater operator comfort ✓ Blue leather cover | One Size | <div>.5mm Pb</div> |
| | Length: 35cm | |



Disposable Lead Gloves

The gloves' thinner texture provides enhanced tactile sensitivity and attenuates up to 50% more direct beam radiation than other leading products. It is powder, latex and lead-free, while the anatomically correct curved fingers reduce hand fatigue.



| PRODUCT OVERVIEW | PROTECTION FOR: | SIZE OPTIONS |
|---|--|--|
| IDEAL FOR HIGH EXPOSURE PROCEDURES <ul style="list-style-type: none"> ✓ Ablations ✓ PTCAs ✓ Biliary drainage ✓ ICDs ✓ Stents ✓ TIPS ✓ Pacemakers ✓ Vascular procedures ✓ Lead extractions | <ul style="list-style-type: none"> • Radiology • Cath/EP Lab • Orthopedics • Endoscopy • Urological procedures • Special procedures • O.R. • Pain management | <div>6 - 9</div> <div>including half sizes</div> |
| | | THICKNESS |
| | | <div>.22mm</div> <div>.30mm</div> |