

Table Cover 1

Basic lower body shield style. This versatile design fits any operating tables including those smaller and narrow sizes.

PRODUCT OVERVIEW

- ✓ Velcro closure system
- ✓ Protective section height/ length (lead panels all together as full protection shield): 72 cm
- ✓ Full protective shield dimensions: 65cm x 72cm
- ✓ Standard protection level supplied: .5mm Pb
- ✓ Lead panel overlapping: 2.5cm
- ✓ Individual panel width: 15cm
- ✓ Main section: 45cmx65cm (other dimensions available on request)
- ✓ Total weight: about 5 kg each – for both of the 2 separate sections
- ✓ Only available in LIGHT lead



Table Cover 2

PRODUCT OVERVIEW

- ✓ Operating table bottom positioning
- ✓ Easy positioning in between the protective lateral side shields with Velcro attachments
- ✓ Only available in LIGHT lead
- ✓ Standard protection level supplied: .5mm Pb
- ✓ Protection shield section dimensions: 52.5x72 cm
- ✓ Panel overlapping: 2.5cm
- ✓ Individual panel width: 15cm
- ✓ Panel height: 72cm
- ✓ Velcro on both sides of the main/ holding section
- ✓ Main/ holding section: 52.5cmx45cm
- ✓ Total weight: about 3 kg roughly

