

Hyposcopic Apron

The table protection screens reduce operator exposure to radiation emissions from below the table and are supplied with a removable upper screen to protect against scatter radiation from the patient's body. Made from a plastic coated lead rubber, the screen is a flexible solution to operator protection.



PRODUCT OVERVIEW

- ✓ Flexible lower screen of plastic coated lead rubber - 0.5mm Pb
- ✓ Mounted on easy remove arm
- ✓ Additional rubber strips can be added to suit alternative size requirements
- ✓ Removable upper screen to protect against scatter radiation from the patient's body
- ✓ Dimensions: Lower screen: 63x75cm (WxH)
Upper screen: 63x22.5cm(WxH)

