

Single Front Apron

(MODEL 07)



3
YEAR
GUARANTEE



Customised
Just For YOU

The Single Front Apron is suitable where protection is required to the front only. While providing complete protection for the user, there is also a choice of three closure options depending on user preference. This design offers complete flexibility and ease of movement, thus ensuring the highest level of comfort at all times.



Velcro closure – Side



Closure A: Velcro flap between the shoulders



Closure B: Buckle with fully adjustable strap



Closure C: Elasticated Velcro fastening provides individual fit

PRODUCT OVERVIEW

- ✓ Greater air circulation
- ✓ Padded shoulders
- ✓ Soft outer material to allow for complete flexibility
- ✓ Even transfer of weight between shoulders and hips to reduce fatigue
- ✓ Made to individual measurements
- ✓ Antibacterial surface

AmRay OPTIONS

- ✓ Monogram
- ✓ Logo on pocket
- ✓ Thyroid collar clipped on
- ✓ Single Front Fluoroscopy Apron available on request



MATERIALS

Light Lead

Light Plus

Supreme Lead-Free


LEAD EQUIVALENCE

.25mm Pb

.35mm Pb

.5mm Pb

PROTECTION



Front protection

SIZE OPTIONS

S

M

L

XL

Other sizes available on request
Please refer to size chart p. 12