

# Semi-Wrap Apron

(MODEL 01)



The Semi-Wrap Apron is a one-piece design affording the user maximum protection and comfort.

Ideally suited for extended procedures such as cardiology, neurology and urology, the 50/25 apron is praised among medical professionals for the significant levels of comfort they enjoy when wearing this ergonomically designed purpose-made apron.



Velcro closure

## PRODUCT OVERVIEW

- ✓ Greater air circulation
- ✓ Padded shoulders
- ✓ Soft outer material to allow for complete flexibility
- ✓ Even transfer of weight between shoulders and hips to reduce fatigue
- ✓ Made to individual measurements
- ✓ Antibacterial surface

## AmRay OPTIONS

- ✓ Monogram
- ✓ Logo on pocket
- ✓ Thyroid collar clipped on
- ✓ Support belt
- ✓ Buckle or Velcro closure



## MATERIALS

Light Lead

Light Plus

Supreme Lead-Free

## LEAD EQUIVALENCE

.35/.25 .35/.35

.50/.25 .50/.35

.70/.25 .70/.35

## PROTECTION



Full protection to front, back and sides with 35cm overlap

## SIZE OPTIONS

XS

S

M

L

XL

2XL

Other sizes available on request  
Please refer to size chart p. 12