

Drop-Off Theatre Apron

(MODEL 08)



The Drop-Off Theatre apron is a single front apron with an instant release Velcro fastening to allow quick and easy removal without breaking the sterile field.



PRODUCT OVERVIEW
<ul style="list-style-type: none"> ✓ Greater air circulation ✓ Padded shoulders ✓ Soft outer material to allow for complete flexibility ✓ Even transfer of weight between shoulders and hips to reduce fatigue ✓ Made to individual measurements ✓ Antibacterial surface

AmRay OPTIONS
<ul style="list-style-type: none"> Monogram Thyroid collar clipped on Velcro closure For Double-Sided Drop-Off Apron, see page 26

MATERIALS
Light Lead
Light Plus
Supreme Lead-Free

LEAD EQUIVALENCE
.25mm Pb
.35mm Pb
.5mm Pb

PROTECTION
Front protection

SIZE OPTIONS	
S	M
L	XL
Other sizes available on request Please refer to size chart p. 12	