













SK THERMO, round Ø 5 mm

Product picture	1:1	Product number	Description	Quality	Recommended speed/rpm	
		A21080052080	SK THERMO 5 mm round	coarse	7.000–10.000 (max. 11.000)	5 x 10*
		A21080052150	SK THERMO 5 mm round	medium	7.000–10.000 (max. 11.000)	5 x 10
		A21080052320	SK THERMO 5 mm round	fine	7.000–10.000 (max. 11.000)	5 x 10





SK THERMO, round Ø 7 mm

Product picture	1:1	Product number	Description	Quality	Recommended speed/rpm	
		A21080072080	SK THERMO 7 mm round	coarse	7.000–10.000 (max. 11.000)	5 x 10
		A21080072150	SK THERMO 7 mm round	medium	7.000–10.000 (max. 11.000)	5 x 10
		A21080072320	SK THERMO 7 mm round	fine	7.000–10.000 (max. 11.000)	5 x 10

SK THERMO, round Ø 10 mm





Product picture	1:1	Product number	Description	Quality	Recommended speed/rpm	
		A21080102080	SK THERMO 10 mm round	coarse	7.000–10.000 (max. 11.000)	5 x 10
		A21080102150	SK THERMO 10 mm round	medium	7.000–10.000 (max. 11.000)	5 x 10
		A21080102320	SK THERMO 10 mm round	fine	7.000–10.000 (max. 11.000)	5 x 10

SK THERMO, round Ø 13 mm

Product picture	1:1	Product number	Description	Quality	Recommended speed/rpm	
		A21080132080	SK THERMO 13 mm round	coarse	7.000–10.000 (max. 11.000)	5 x 10
		A21080132150	SK THERMO 13 mm round	medium	7.000–10.000 (max. 11.000)	5 x 10
		A21080132320	SK THERMO 13 mm round	fine	7.000–10.000 (max. 11.000)	5 x 10

* Every box contains 5 PE bags with 10 abrasive caps each.

SK THERMO, round Ø 16 mm

Product picture	1:1	Product number	Description	Quality	empf. Drehzahl U/min	
		A21080162080	SK THERMO 16 mm round	coarse	7.000 – 10.000 (max. 11.000)	5 x 10
		A21080162150	SK THERMO 16 mm round	medium	7.000 – 10.000 (max. 11.000)	5 x 10
		A21080162320	SK THERMO 16 mm round	fine	7.000 – 10.000 (max. 11.000)	5 x 10

Precise abrasive cap carriers for an exact fit

The SK THERMO fits precisely on the LUKAS abrasive cap carriers, allowing it to unfold its full performance. During the treatment, the cap sits tightly on the carrier. After use, it can easily be exchanged.



See page 13 for matching abrasive cap carriers