

Classic Apron

(MODEL 09)



The **AmRay Classic** is a low budget, functional single front apron. It is designed for high usage environment and withstands the rigours of everyday use. Ideal for short periods of use only.



PRODUCT OVERVIEW

- ✓ Greater air circulation
- ✓ Even transfer of weight between shoulders and hips to reduce fatigue
- ✓ Made to individual measurements
- ✓ Low-cost functional apron
- ✓ Buckle closure

AmRay OPTIONS


- ✓ Monogram
- ✓ **A-Track**
AmRay Inventory Management

Colour OPTIONS

Only available in nylon
Only available in 2 colours



Blue Nylon



Navy Nylon

MATERIALS

Light Lead


LEAD EQUIVALENCE

.25mm Pb

.35mm Pb

.5mm Pb

PROTECTION



Front protection

SIZE OPTIONS

S

M

L

XL

Other sizes available on request
Please refer to size chart p. 12