

Backsaver

(MODEL 02)

The Backsaver offers all-round protection with less weight and is designed with adjustable elasticated back panels for back support and individual comfort.



Velcro closure



PRODUCT OVERVIEW

- ✓ Greater air circulation
- ✓ Padded shoulders
- ✓ Soft outer material to allow for complete flexibility
- ✓ Even transfer of weight between shoulders and hips to reduce fatigue
- ✓ Made to individual measurements
- ✓ Antibacterial surface

AmRay OPTIONS

- ✓ Monogram
- ✓ Logo on pocket
- ✓ Thyroid collar clipped on
- ✓ Velcro closure
- ✓  A-Track
AmRay Inventory Management

MATERIALS

Light Lead

Light Plus

Supreme Lead-Free

LEAD EQUIVALENCE

.25/.25

.35/.25

.35/.35

.50/.25

.50/.35

PROTECTION



Front, back and side protection

SIZE OPTIONS

XS

S

M

L

XL

2XL

Other sizes available on request
Please refer to size chart p. 12