

Mobile Hangers

As you have invested in a personal protective apron, it is important that it is maintained properly so that the maximum lifespan is achieved. With the versatility of a mobile hanger, all sets of protective garments can be easily stored in the proper way. All **AmRay** mobile hangers are designed to be easily moved with or without hung garments.



Mobile Multi-hanger (MODEL AMMH6)

PRODUCT OVERVIEW	DIMENSIONS
<ul style="list-style-type: none"> ✓ 6 aprons mobile hanger ✓ Revolving supports ✓ Lockable castors for ease of movement ✓ Made in powder-coated steel 	<p>Height: 144cm</p> <p>Width: 48cm</p> <p>Length of arm: 57cm</p>



Mobile hanger (MODEL AMMH10)

PRODUCT OVERVIEW	DIMENSIONS
<ul style="list-style-type: none"> ✓ 10 aprons mobile hanger ✓ Revolving supports ✓ Lockable castors for ease of movement ✓ Made in powder-coated steel 	<p>Height: 132cm</p> <p>Width: 81cm</p> <p>Length: 58cm</p>



Mobile Double-hanger (MODEL AMMH1010)

PRODUCT OVERVIEW	DIMENSIONS
<ul style="list-style-type: none"> ✓ 10 aprons mobile hanger ✓ Revolving supports ✓ Lockable castors for ease of movement ✓ Made in powder-coated steel 	<p>Height: 135cm</p> <p>Width: 99.6cm</p> <p>Depth: 64.5cm</p> <p>Arm length: 50cm</p>



Mobile Multi-hanger (MODEL AMMH55)

PRODUCT OVERVIEW	DIMENSIONS
<ul style="list-style-type: none"> ✓ Stores up to 5 lead aprons, which can be a combination of two-piece aprons, like vest and skirt, or one-piece aprons. ✓ Lockable castors for ease of movement ✓ Made in powder-coated steel 	<p>Height: 135cm</p> <p>Width: 52.3cm</p> <p>Depth: 64.5cm</p> <p>Arm length: 50cm</p>