

AV4036 - AV4022

Visit couch 2 sections 1 joint, fixed height

Givas visit couches models AV4036 and AV4022 are characterised by high level of convenience and for their comfortable dimensions. The materials that have been carefully chosen, together with the linear structure help operators during cleaning and maintenance operations. Upper structure made of painted steel while the legs are made of chromed steel and connected to the frame through meltings of alloy. The gas-spring for backrest adjustment is placed on the upper structure. Lying parts composed by a support of wood conglomerate, padded with a waterproof polyurethane slab and upholstered with fireproof leatherette fabric, class 11M. The controls for backrest are placed on both sides of the couch: gas-spring for AV4036 or double rack (6 positions) for AV4022. Foot section is fixed. The upholstery is made of plastic washable material, certificate of waterproof class 11M, seamless and easily disinfected.



AV4036
Backrest adjustable through gas-spring with control on both sides of the bed

AV4022
Backrest adjustable through a double semiautomatic rack (6 positions)



Available accessories



AV4097B
Paper support made of chromed steel, applied to bed legs, complete with clamps made of plastic



AV4088
Couple of lateral and universal profiles, for accessories coupling



AV4098
I.V stand with two hooks made of chromed steel, complete with clamp, adjustable in height and applicable to the lateral rails

Dati tecnici Technical features

DATI DIMENSIONALI / DIMENSIONS	
Dimensioni totali / Total dimensions	L 190 / P 56 / H 80 cm
Dimensioni piano di coricamento / Lying area dimensions	L 190 / P 56 cm
Dimensioni schienale / Back rest dimensions	L 56 / P 56 cm
Dimensione gambale / Foot rest dimensions	L 134 / H 56 cm
Inclinazione schienale / Back rest tilt	70°
Portata / Safe operating load	150 kg
Peso / Weight	29 Kg
Dispositivo medico marcato CE DIR.93/42 appartenente alla classe 1 / Medical device Class 1 following dir. 93/42 CEE	
Prog. di sistema attribuito al D.M.66589 Classif. CND Y181207 / Registration no. at MOH 66589 / CND Y18120	