

AT4054 - AT4056

Visit couch 2 sections 1 joint, fixed height



The series of therapy couches Cardiomass Givas AT4054 and AT4056, offers professional couches with high versatility, suitable for all types of medical and estetical treatment, therapy or massage. The simple adjustment of back section through crank, facilitate medical staff for all eventualities, assuring comfort and safety for the patient during the treatment.


The structure of the couch is made of varnished solid wood , robust and resistant to moisture.

The back section can be adjusted through a rack. The lying surface is padded with soft polyurethane; the external covering is made of fireproof plastic material with class 1IM and easily disinfected. Seams are double to be more safe and tear resistant.




AT4054

AT4054



AT4056



Available accessories



AT4080

Element with face hole for therapy bed



AT4097C

Paper support for cardiomass bed



AT4086

Wedge-shaped cushion

Dati tecnici Technical features

DATI DIMENSIONALI / DIMENSIONS		
AT4054	Dimensioni totali / Total dimensions	L 190 / P 76 / H 80 cm
AT4056	Dimensioni totali / Total dimensions	L 190 / P 76 / H 66 cm
Dimensioni piano di coricamento / Lying area dimensions		L 190 / P 76 cm
Dimensioni schienale / Back rest dimensions		L 60 / P 76 cm
Dimensione gambale / Foot rest dimensions		L 130 / H 76 cm
Inclinazione schienale / Back rest tilt		32°
Portata / Capacity		150 kg
Peso / Weight		42 Kg
Dispositivo medico marcato CE DIR.93/42 appartenente alla classe 1 / Medical device Class 1 following dir. 93/42 CEE		
Prog. di sistema attribuito al D.M.66589 Classif. CND Y181207 / Registration no. at MOH 66589 / CND Y18120		