

## **CHILD COT**

### **MODEL DC2000-C**



#### **Product Features**

- Steel base with ventilation holes
- Available with wooden head/head board with cartoon picture
- Tubular steel pole head & leg frame
- Scratch and chip resistant epoxy powder coating
- Safety side rails of filling not exceeding 102mm spacing with rubber bushes along 2 chrome rods fitted on leg frame eliminating friction noises
- Auto engage side rail
- Mobile on 4 x 100mm swivel castors with diagonally fitted brakes

#### **Technical Specifications**

- Approx. Dim (mm): 1400(L) x 850(W) x 665(H) (from floor to mattress platform)
- Mattress platform to top of side rail: 590mm (H)
  - \*Optional accessories: Mediflex mattress
- Mattress Approx. Dim (mm): 1285(L) x 762(W) x 102(Thick)
  - \*Available for adjustable backrest by mechanical crank winding mechanism in model **DC2010-C** (Single fowler)