

## ✚✚✚ PVC Anesthesia Mask pvc

1. Made of medical grade PVC, clear and soft.
2. Seal tight against mouth.
3. Color - coded to easily identify sizes.

Item No.	Size (#)
AMP0	0
AMP1	1
AMP2	2
AMP3	3
AMP4	4
AMP5	5
AMP6	6



## ✚✚✚ Silicone Manual Resuscitator

1. The product includes silicone resuscitator (1), silicone face mask (1), oxygen tube (1), reservoir bag (1), airways (3) and mouth opener (1). Mainly used for pulmonary resuscitation and artificial respiration.
2. Three sizes available, suitable for adult, child and infant.
3. Can be repeated use, 100% latex free.

Item No.	Size
SMR-L	Adult
SMR-M	Pediatric
SMR-S	Infant



## ✚✚✚ PVC Manual Resuscitator PVC


1. This product is intended for pulmonary resuscitation for patients.
2. The accessories includes PVC mask (1), oxygen tubing (1), reservoir bag (1), opener (1) and airways (3).
3. For single use only, 100% latex free.

Item No.	Size
PMR-L	Adult
PMR-M	Pediatric
PMR-S	Infant



## ✚✚✚ Guedel Airway

1. Made of semi-rigid, non-toxic polyethylene.
2. Guedel style enclosed center channel.
3. Color coded to easily identify sizes.



Item No.	Size (#)	Length (mm)	Color Code
GA40	000	40	Pink
GA50	00	50	Blue
GA60	0	60	Black
GA70	1	70	White
GA80	2	80	Green
GA90	3	90	Yellow
GA100	4	100	Red
GA110	5	110	Light Blue
GA120	6	120	Orange

